

## **MLA1010B Disposable ECG Electrodes (1000)**

### *Accessories*

A pack of 1000 disposable Ag/AgCl ECG electrodes for use with snap connect lead wires MLA0313, MLA0315, MLA2503 and MLA2505. The electrodes can be used for ECG, EOG and EMG.

Preparation: Before attaching the electrode, gently abrade the skin using abrasive gel (MLA1093). This reduces electrical resistance and can improve the quality of recordings.



#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA1010B Disposable ECG Electrodes (1000)

For use with:

MLA0313 Lead Wires (3-snap on)

MLA0315 Lead Wires (5-snap on)

MLA2503 Shielded Lead Wires (3-snap on)

MLA2505 Shielded Lead Wires (5-snap on)