

MLA1093 Abrasive Gel

Accessories

The MLA1093 Abrasive Gel 4oz tube can be used to prepare an EEG, EOG or EMG recording site. The mild abrasive in the gel helps lower impedances at the electrode site with minimal discomfort and irritation. Apply gel to electrode site with gauze pad or cotton swab..



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1093 Abrasive Gel

A pack of 3 tubes of Abrasive Gel is also available: MLA1093B