

MLA1093B Abrasive Gel (3 Pack)

Accessories

The Abrasive Gel (3 x 4 oz tubes) can be used to prepare an EEG, ECG or EMG recording site. The mild abrasive in the gel helps lower impedances at the electrode site with minimal discomfort and irritation. Apply gel to electrode site with gauze pad or cotton swab, rubbing gently into skin.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1093B Abrasive Gel (3 pack)