

## **MLA255 Differential Pod Input Cable**

### *Accessories*

Cable with an 8 pin DIN connector to three alligator clips suitable as a recording cable in experiments with the MLT012 Nerve Chamber. Provides differential input to PowerLab models with Pod ports or the FE305 Pod Expander.



#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA255 Differential Pod Input Cable

For use with:

MLT012 Nerve Chamber

FE305 Pod Expander

Any PowerLab with Pod ports