

## **MLA304 Disposable Filters**

### *Accessories*

---

A Pack of 50 Disposable Filters protect against cross contamination without impairing airflow. MLA1026 Reusable Mouthpieces and MLA1011A Clean Bore Tubing & Flow Head Adapter can be used with the filters. Intended for single human use only — do not wash or sterilize.



Dimensions:

Subject end 31 mm (1.2") outer diameter

Tubing end 34 mm (1.3") outer diameter

#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA304 Disposable Filters

For use with:

MLA1026 Reusable Mouthpieces

MLA1011A Clean Bore Tubing and Flow Head Adapter