

## **MLA5530 Calibration Syringe 3 L**

### *Accessories*

The MLA5530 3 L Calibration Syringe provides an easy and reliable method of calibrating respiratory equipment. Other volumes ranging from 10 mL to 7 L are available on request.



When calibrating, a steady rate of movement of the plunger should be maintained. High rates of movement and sudden jolts may cause unwanted interference on the recorded signal.

#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA5530 Calibration Syringe 3 L

For use with:

MLT1000L Flow Head