

## **MLT117F GSR Finger Electrodes (MR Safe)**



### *Transducer Series*

MRI safe, bipolar finger electrodes for measuring galvanic skin response. The silver/silver chloride contoured finger electrodes are supplied with Velcro straps, a 5 m cable (16') and 2 x 4 mm Shrouded connectors for connection to the FE116 GSR Amp.



It is recommended that these electrodes are used with MLA1095 Electrode Paste, to improve electrode to skin contact and GSR recordings.

#### **Caution**

Read "Statement of Intended Use" on our website or in "Getting Started with PowerLab" before use.

#### **Ordering Information:**

MLT117F GSR Finger Electrodes

For use with:  
FE116 GSR Amp